


Improved Efficiency of Medical Imaging Probe



S. L. Newton, I. Styles, M. M. Britton, A. F. A. Peacock PSIBS, School of Chemistry, University of Birmingham, B17 2TT SLN286@bham.ac.uk

The Peacock Group has developed a new MRI imaging probe with the potential to drastically improve the identification of disease compared to the current standard. This new design uses a miniature protein which binds to an active metal.


Medical Standard

MRI is a technique commonly used in hospitals to look inside the body without the need for invasive surgery.

Contrast agents can be injected to allow for easier diagnosis of disease.


Factors affecting contrast agents:

- 1. Rate of tumbling (peptide length)
- 2. Association with water




Peptide design

Each amino acid was selected to drive the formation of the coiled coil structure or to form a metal binding site for an MRI active metal. In this system the metal is gadolinium.



Peptide strands Ac-G-IAANEWKDAAIEQKIAAIEQKIAAIEQKIAAIEQK-G-NH2

1. Investigating peptide length



An increase in peptide length causes an increase in folding. This means the peptide is more stable and can bind a metal more tightly. This is important as on its own gadolinium is highly toxic.


MRI has shown that the efficiency increases peptide length as However increases. there is no difference between the medium and long peptide

Peptide	Efficiency r2 / mM-1 s-1	
Medical Standard	16.8 ± 0.5	
Short	54.1 ± 1.0	
Medium	92.5 ± 11.0	
Long	82.7 ± 1.9	

2. Investigating water association

Water can associate with the outside of the peptide or with the metal in the middle of the peptide. Here we investigate what happens if more water bonding groups are added to the outside of the structure.

The modification in structure causes no significant change in peptide folding or how tightly the metal is bound.



But the new peptide has more water coordinated to the metal centre.

Computational simulations show one water molecule in the centre of the coil

and many coordinated to the outside.

This results in an increase in efficiency of the new peptide.

	1, 21/2	•
•		*
		. 6
	10	

bind water



Longer length = slower tumbling = higher efficiency

More H-bonding groups = more water binding = higher efficiency

Acknowledgments

would like to thank my supervisors for help planning the experiments and instrument training, John Wilke for help with MD work, thanks to AWM for use of equipment, the School of Chemistry and PSIBS and the EPSRC for funding.

References

P. Caravan. Acc. Chem. Res. 42, 851 (2009). M.R. Berwick et al. J. Am. Chem. Soc. 136, 1166 (2014).

D. K. Wilkins et al. Biochem. 38, 16424 (1999).

A.D. Sherry et al, Current Opinions in Chemical Biology, 2013, 17: 167-174



BIRMINGHAM









